

Best[™] nomos[®]

Compact Sonalis[®] Ultrasound System



The Best[™] Sonalis[®] Ultrasound Imaging System provides superior visualization of HDR, LDR, RF or Cryosurgical procedures. Our patented SimulView[™] Technology provides simultaneous “live” views of the prostate in both planes, thereby increasing treatment accuracy and precision.

Additional Features

- Capable of supporting over 20 different probe configurations, this system can be used in almost every discipline
- 21.5" touch screen technology
- Superior HD image resolution for improved procedure accuracy, speed and physician confidence level
- Patented probe design
- Sagittal array provides for 140 mm length of view encompassing the bladder, prostate and perineal tissue
- Advanced modular software design provides for future upgrade path
- Convertible to desktop or portable without losing any functionality
- Dual battery backup built-in to protect from power loss



Best Nomos, 1 Best Drive, Pittsburgh, PA 15202 USA

tel: 412 312 6700 800 70 NOMOS fax: 412 312 6701 www.nomos.com

*Certain products shown are not available for sale currently. TeamBest Global Companies ©2024