



# Best™ *medical international*

A TEAMBEST GLOBAL COMPANY

TEAMBEST GLOBAL®



www.teambest.com



## Best GammaBeam™ 100/300 Equinox™ Teletherapy System with Avanza™ 6D Patient Positioning Table

With **NEW** Multi-Leaf Collimator for 80 and 100 cm SAD units—IMRT, IGRT, SRS, SBRT and Tomotherapy capable with ActiveRx



Best™ Gammacell® 1000/3000 Blood & Research Irradiator

## UPGRADE KIT for all old Theratron units, 80 or 100 cm including IMRT capabilities w/built in or external MLC



### UPGRADE includes:

- Removing all old controls, electronics and installing a new control system and covers
- Replacing the old collimator system with the new Equinox collimator
- Replacing the old treatment table with the new Avanza™ Table
- Retaining the head rotation capability is optional

### UPGRADE features:

- Calculated Arc Speed
- Graphical Control System
- Asymmetric Jaws (optional)
- R&V System Ready (optional)
- Service Log Files
- On-Board Verification
- Motorized Wedge (optional)
- Collision Detection (optional)



Best™ Raycell Mk2 Blood & Research Irradiator



Best™ Raycell Mk1 Blood & Research Irradiator

Best Theratronics Ltd., 413 March Road, Ottawa, Ontario K2K 0E4 Canada  
tel: 613 591 2100 866 792 8598 www.theratronics.com www.teambest.com

\*Certain products shown are not available for sale currently. TeamBest Global Companies ©2024